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EN Treatment of discolored pulpless teeth

FORMULA

100 % carbamide peroxide.

INDICATIONS

Treatment of discolored dentin in pulpless teeth.

CONTRAINDICATIONS

External application on the teeth is formally contraindicated.

PROPERTIES

When in contact with physiological fluids, carbamide peroxide releases oxygen. This occurs in an almost neutral medium.

Its main properties are due to the oxidizing power of active oxygen which is: a bleaching agent, a bactericide and a detergent by gas evolution.

INSTRUCTIONS FOR USE

Before starting treatment, make sure that all canals have been thoroughly filled with a hard paste.

- Clear the pulp chamber to the floor
- Make a paste by mixing two spoonfuls of powder with 1 drop of glycerol to a firm consistency
- Carefully cover the floor and walls of the cavity with the obtained paste
- Cover with a cotton pellet and seal with cement

Re-examine the patient after 4 or 5 days. If the shade of the tooth is light enough, seal. If not, repeat treatment.

It is best to spread the treatment over several sessions, so as not to exceed the desired shade.

SPECIAL WARNINGS AND PRECAUTIONS FOR USE

C - Corrosive

Contains carbamide peroxide.

Causes burns.

In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

Wear suitable protective clothing, gloves and eye/face protection.

In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

STORAGE

Keep away from light in its original packaging and do not store above 25°C.

Close the container tightly after use.

PRESENTATION

5 g bottle + 1 measuring spoon.

Kit containing :

- 1 x 5 g bottle of powder
- 1 x 5 ml flask of glycerol
- 1 measuring spoon

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