F

COMPOSITION

Eucalyptol

PROPRIÉTÉS

Solvant de la gutta-percha.

INDICATIONS

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Solvent for gutta percha.

INDICATIONS

Fitting of the primary gutta percha point in root-canal filling techniques using the lateral condensation method.
Retreatment after root canal filling by gutta percha condensation.

INSTRUCTIONS FOR USE

- FITTING OF THE PRIMARY POINT (LATERAL CONDENSATION)
With the point-fitting technique, the solvent softens the surface of the primary point so that it adapts perfectly to the irregularities and contour of the root canal.
First, the dentist selects the primary point depending on the type and size of the root canal: the apical tip of the point is dipped into the solvent.
Second, the point is inserted into the canal over the whole working length using a slight up and down motion.
The point is then withdrawn from the canal and coated with cement sealer.
It is reinserted by gentle pumping movements so as to carry the cement to the apex.
Then, using a solvent-lubricated spreader, the point is condensed laterally.
- RETREATMENT:
Place two drops of solvent into the pulp chamber.
Then, using a manual root canal instrument, remove the softened gutta percha until the instrument reaches the apical foramen.
Add some fresh solvent from time to time and rinse alternately with sodium hypochlorite.

SPECIAL WARNINGS AND PRECAUTIONS FOR USE

Flammable.

Avoid contact with skin and eyes.

Do not swallow.

STORAGE

Keep at room temperature in the original container, away from light and moisture.

PRESENTATION

13 ml bottle.

For professional use only in the practice of dentistry.

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Guttasolv
Solvent for gutta percha

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