FR

PROPRIÉTÉS

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EN

PROPERTIES

PROVISCELL, an eugenol-free temporary cement is simple to apply. PROVISCELL has no influence on the hardening of resin-based sealing or filling materials.

This product is in accordance to the standard ISO 3107.

INDICATION

Temporary sealing.

CONTRAINDICATION

Do not use PROVISCELL if the patient is know to be allergic to fatty acid and colophany.

DIRECTIONS FOR USE

• Tubes

PROVISCELL should be used at room temperature (18°C- 23°C). In case of lower temperatures (a cool glass slab is used), the working time will be longer.

- Squeeze equal amounts of base and catalyst into a mixing pad or glass plate.
- Mix the two pastes for approximately 20 seconds to a homogeneous consistency. Use a plastic or rustproof metal spatula for mixing. To obtain a softer material and a longer working time add a small amount of petrolatum.
- A harder and faster setting material is obtained by adding zinc oxide.
- Apply a thin layer of PROVISCELL inside the temporary crown or bridge and cement it to the prepared tooth.
- PROVISCELL (tube) remains workable for about 1 1/2-2 minutes and sets in the mouth within 4-6 minutes.
- 3-5 minutes after PROVISCELL has set, remove excess material.

• Syringes

- Remove the cap by turning 90° anticlockwise and discard it.
- Attach a mixing tip to the syringe. Line up the tip notch and the syringe notch.
- Turn the mixing tip 90° clockwise until it locks.
- The cement is automatically mixed in the tip and can be applied directly.
- At first use, extrude pea-sized amount of material and discard it.
- Apply the mixed cement into the dried temporary restoration.
- Fix the temporary restoration onto the dried prepared core, exerting gentle pressure.
- After setting (approximately 4 minutes), remove excess cement with a probe.

Leave the used mixing tip on the syringe to serve as a cap until next use.

SPECIAL WARNINGS

In case of exposed pulp or very low dentin thickness, it is recommended to use a pulp cap before.

N - Dangerous for the environment

Very toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

Xi - Irritating to eyes, respiratory system and skin.

May cause sensitisation by skin contact.

Avoid contact with skin and eyes.

In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

Wear suitable gloves and eye/face protection.

STORAGE

Keep the tubes at temperatures between 4°C and 23°C and protect from moisture.

Keep the syringes at temperatures between 15°C and 25°C. Do not keep in refrigerator. Protect from direct sunlight.

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