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Bioactive mineral root canal sealer

COMPOSITION

Powder based on tricalcium silicate, zirconium oxide and excipients. Aqueous solution of calcium chloride and excipients.

PROPERTIES

BioRoot™ RCS is a bioactive mineral root canal sealer based on Septodont patented and innovative mineral micro-aggregate chemistry "Active Biosilicate Technology" that offers:

- 1) Biocompatibility: high mineral purity and monomer free formulation reducing the risk of adverse tissue reaction.
- 2) Bioactive properties: hydroxyapatite formation at the tooth-sealer interface and mineralization of dentinal structure.
- 3) Alkaline pH.
- 4) Sealing properties: crystallization of the material inside the dentin tubules creating a tight seal.
- 5) Retreatable: if a retreatment is needed, BioRoot™ RCS can be easily removed from the root canal.

INDICATIONS

Permanent root canal filling in combination with gutta-percha points in vital or necrotic pulp cases or after a retreatment procedure.

BioRoot™ RCS is suitable for use in single cone technique or cold lateral condensation.

CONTRAINDICTION & PRECAUTION

Allergy to one of the components.

ADVERSE EFFECTS

No adverse effects known.

MIXING INSTRUCTIONS

The mix must be made extemporaneously.

- 1/ With the spoon supplied in the box, select the spoon cavity according to the expected product quantity.
- 2/ Collect the powder. Put a level spoonful of the powder on the mixing pad.
- 3/ Detach a mixing solution single-dose container. Twist cap to open.
- 4/ Pour the following number of drops on the mixing pad:
 - 3 drops for the small spoon cavity
 - 5 drops for the large spoon cavity
- 5/ Prepare the root canal sealer by progressively adding powder to the liquid. Mix until obtaining a smooth paste (about 60 seconds).
- 6/ Immediately rinse and clean the instruments to remove any residual material.
- 7/ BioRoot™ RCS has a minimum working time of 10 minutes and a maximum setting time of 4 hours.

PLACEMENT INTO THE ROOT CANAL

- 1/ Perform the root canal cleaning and shaping procedure using standard endodontic procedures.
- 2/ Select standardized gutta-percha master cone and check that it fits snugly at the working length.
- 3/ Dry the canal with paper points.
- 4/ Prepare BioRoot™ RCS according to the above mentioned mixing instructions.
- 5/ Apply a coating of BioRoot™ RCS onto the canal walls using either a paper point or the gutta-percha point.
- 6/ Complete the obturation by inserting the gutta-percha master cone previously coated with BioRoot™ RCS (single cone technique) or several coated gutta-percha points (lateral condensation technique).
- 7/ Assess the quality of the root canal filling with a radiograph.
- 8/ Use a hot instrument to sever the protruding gutta-percha.
- 9/ Perform the temporary restoration using any temporary filling material.
- 10/ If a post placement is required, it should be performed after complete hardening of BioRoot™ RCS, i.e. between 1 week and 1 month.

A radiograph should be taken to ascertain the persistence of 3 to 5 mm of the apical root filling.

BioRoot™ RCS REMOVAL

Since BioRoot™ RCS is used in combination with gutta-percha points, it can be removed using conventional removal techniques.

WARNINGS & PRECAUTIONS FOR USE

- Ensure that the rubber dam is properly placed so as to completely isolate the operating field.
- In case of a wide apical foramen, an apexification of the root should precede the root canal filling.

CONSERVATION

Store in a dry place.

Single-dose containers must be used immediately after opening. Any remaining material should be discarded.

Use the powder within 6 months after first opening.

PRESENTATIONS

BioRoot™ RCS:

- 1 bottle of 15 g powder
- 35 liquid single-dose containers
- 1 spoon

BioRoot™ RCS Liquid single-dose Container:

- 20 liquid single-dose containers

Made by:

SEPTODONT

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94100 Saint-Maur-des-Fossés - France

Tel : 33 (1) 49 76 70 00

Distributed in USA by:

SEPTODONT

Louisville CO 80027 - USA

**For Professional Dental Use Only
Rx Only**

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