

SCALPEL-DESIGNED BEVEL



EASIER
INSERTION



REDUCED
DEFLECTION

SEPTOJECT EVOLUTION

GREATER CONTROL, LESS DISCOMFORT.



SEPTOJECT EVOLUTION

THE INNOVATIVE NEEDLE WITH A PATENTED SCALPEL-DESIGNED BEVEL.

Septoject Evolution is engineered for a smoother penetration with less tissue displacement. This means less discomfort for your patients, even when used for multiple injections. For you, the need for less force, combined with significantly reduced deflection, brings you better control and accuracy.

SCALPEL-DESIGNED BEVEL



EASIER
INSERTION



REDUCED
DEFLECTION



GREATER CONTROL REDUCED DISCOMFORT


• Easier Insertion

• Easy Penetration:

29% less force is required for insertion using Septoject Evolution.

• Ideal for multiple injections:

The fourth insertion with Septoject Evolution still requires less force than the first insertion with a standard needle.





• Reduced Deflection

• Better accuracy:


Septoject Evolution produces more limited deflection.

Comparative Deflection (mm) of Needle After Penetration Through a 10mm Block of Silicone²



Average Distance to Target (mm).
The Closer to 0, the Less Deflection

DESIGNED TO MAKE IT EASY FOR YOU



INFILTRATIONS

- Infiltrations (periapical)
- PSA – Posterior Superior Alveolar nerve block
- MSA – Middle Superior Alveolar nerve block
- Palatal Anaesthesia

GAUGE	LENGTH	COLOR
27G	Short 16 mm	Yellow
27G	Short 25 mm	Yellow
30G	Short 16 mm	Blue
30G	Short 25 mm	Blue



INTRALIGAMENTARY INJECTIONS

- Intraligamentary Injections
- Intraseptal Injections
- Intraosseous

GAUGE	LENGTH	COLOR
27G	X-Short 9 mm	Gold
30G	X-Short 9 mm	Purple



If you think a needle is just a needle, think again.
SEPTOJECT EVOLUTION.

